



FootTrak

Motion Analysis Corporation's FootTrak outputs movement analysis of the rearfoot and forefoot and can be used to compare the impact of footwear and orthotics on foot and ankle motion during walking and running.

FootTrak is compatible with all Motion Analysis motion capture systems and requires only basic knowledge of subject set-up and processing.

Features

- Multi-segment foot model
- Rearfoot motion (pronation, supination)
- Ankle motion
- Inversion-eversion
- Pronation velocity
- Simple and intuitive user interface enables detailed, efficient report generation
- Printable and web-ready dashboards



Motion Analysis Corporation
3617 Westwind Boulevard
Santa Rosa, CA USA 95403



707.579.6500
www.motionanalysis.com
sales@motionanalysis.com