




Motion Composer is a software for collating and presenting interactive motion capture Embedded within Cortex, Motion Composer is able to continuously update information data. without exiting to a stand-alone software package, creating a seamless workflow. Users are able to export a user-defined, interactive, graphical report comprising: motion capture data, forceplate data, data, kinetics and kinematics, documents, color video and live web content. Motion Composer EMG reports can be viewed on any computer by using the Motion Viewer presentation player.



Features

- Cortex plugin for presentation authoring
- Presentation tree/property controls for managing presentation resources
- Slide creation tools for laying out presentation slides
- Presentation loading/saving to presentation file format (XML)
- Presentation 'Pack-and-Go' for one-click packaging of viewer and data into a single file
- Freely redistributable viewer for presentation sharing
- Access to all data views available in Cortex
- RTF editor pane for rtf/text create/import/view /edit
- HTML editor pane for HTML/media import/view /edit
- 2D graph pane for graphical information view /analysis

Motion Composer Includes

Motion Composer

Presentation authoring plug-in for Cortex. Author presentations seamlessly from your Cortex session.

Motion Viewer

Presentation player. Package and share presentations easily with this interactive viewer. Suitable for burning to a CD, emailing, or web distribution.

Data Views

Graph, document and multimedia views. Perfect for presentation data.